Swedish Leaf Spring AFO Orthosis Drop Foot Brace

SKU: FTB302

Instructions For Use

Step 1. Prepare the Brace

Unstrap the velcro and place your heel securely in the back of the brace.

Step 2. Adjust the AFO

Mark the length of your foot and trim the AFO with scissors to desired length.

Step 3. Secure Foot

Insert AFO into shoe, then place foot into shoe. Secure Velcro strap around the calf. Trim Velcro strap if needed.



Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do **NOT** dry clean, iron, or bleach.